

Weights and measurements microscopes Labomed

LABOMED
Long, wide, high
Date 7-11-2019



Prima DNT / ENT / GN / OPH	Weight with box	Weight without box	Sizes (L, An, Alt) Without box in cm	sizes (L, An, Alt) with box in cm
Base	70kg	60kg	56 x 54 x 20	70x33x33cm
Column	22kg	10kg	6 x 150cm	164 x 15 x 13 cm
short arm	18kg	14kg	125x8x28cm	158 x 35 x 28 cm
Long arm	23kg	18kg	155x7,5x27cm	176 x 43,5 x 26,5 cm
Short arm with head	22,5kg	16,5kg	110x7,5x27cm	159x42x37cm
Long arm with head	30kg	20,5kg	150x8x28cm	159x42x37cm
Head ENT	7kg	4,6 kg	46x50x24cm	58x43x24cm
Head GN	7kg	4,6 kg	46x50x24cm	58x43x24cm
Head DNT	7kg	4,6 kg	46x50x24cm	58x43x24cm
Head OPH	7kg	4,6 kg	46x50x24cm	58x43x24cm
Head with eyepieces	1kg	0,2kg	4,5 x 5cm	5x5x7cm
Table mounting system	12kg	7kg	65x13x50cm	84x15x13cm

Magna	Weight with box	Weight without box	Sizes (L, W, H) Without box in cm	sizes (L, W, H) with box in cm
Base	91kg	80kg	85x85x20cm	106x42x24cm
Column	47kg	18,5kg	100x7,5x27cm	122x68x45cm
Long arm	47kg	20kg	140x9x28cm	122x68x45cm
Long arm with head	41kg	31kg	167x11x26cm	122x68x45cm
Head	19,4kg	11kg	72x27x40cm	96x61x39cm
Head with eyepieces	1kg	0,2kg	45mm (rond)	5x5x7cm

Weights and measurements microscopes Labomed

Prima Mu	Weight with box	Weight without box	Sizes (L, W, H) Without box in cm	sizes (L, W, H) with box in cm
Base	70	60kg	55x55x20cm	70x33x33cm
Column	13kg	8,5kg	125x7,5x27cm	164x15x13cm
short arm	18kg	14kg	120x8x28cm	159x42x37cm
Short arm with head	23kg	19kg	125x9x29cm	159x42x37cm
Head	7kg	4,6 kg	46x50x24cm	58x43x24cm
Head with eyepieces	1kg	0,2kg	35mm (rond)	5x5x7cm
MuTrainer with mounting system Art. 6138000-255	33kg	25,5kg	85x85x28cm	136x50x38cm
Table mounting system	11kg	6,5kg	65x13x50cm	84x15x13cm

Prima C_CS	Weight with box	Weight without box	Sizes (L, W, H) Without box in cm	sizes (L, W, H) with box in cm
Base	60kg	50kg	56x54x20cm	1x box 67x33x29cm
Column	10kg	8kg	6x150cm	136x51x37cm
Arrm	13kg	10kg	85x85x28cm	136x51x37cm
Head C	6kg	4,7kg	46x25x24cm	53x33x30cm
Head CS	6kg	4,7kg	46x25x24cm	53x33x30cm
Head with eyepieces	1kg	0,2kg	4,5x5cm	5x5x7cm

Weights and measurements microscopes Labomed

Slit Lamps	Weight with box	Weight without box	Width X Depth X Height	
			Sizes (L, An, Alt) Without box in cm	Sizes (L, An, Alt) with box in cm
SLx20 Slit Lamp, 2 step, LED	26kg	20kg	47x42x60cm	58x46x75cm
SLx40 Slit Lamp, 3 step, LED	25kg	19kg	47x42x60cm	58x46x75cm
SLx45 Slit Lamp, 3 step, LED	32kg	24kg	47x42x60cm	73x73x50cm
eVO300 Slit Lamp, 5 step, LED	27kg	21kg	27x36x51cm	58x46x75cm
eVO400 Slit Lamp, 5 step, LED	30kg	23kg	44x39x68cm	73x73x50cm
eVO400 Slit Lamp, 5 step, Halogen	30kg	23kg	44x39x68cm	73x73x50cm
PRIMA eVO500D Slit Lamp, 5 step LED	27kg	21kg	53x50x70cm	70x47x87cm