

GYNECOLOGY



LABOMED® 



BUDGET FAVORITE

PRIMA C

Illumination
27W LED, 75.000 Lux

Vertical movement of arm
100 mm



SUPER ALL-ROUNDER

PRIMA CS

Illumination
27W LED, 75.000 Lux

Vertical movement of arm
± 250 mm



SUPREME CHOICE

PRIMA GN

Illumination
50W LED, 135.000 Lux

Vertical movement of arm
500 mm



OVER-THE-SHOULDER

PRIMA GN OTS

Illumination
50W LED, 135.000 Lux

Vertical movement of arm
500 mm



BUDGET FAVORITE

Compact and versatile colposcope,
high power LED, superb optics
and ergonomic footprint.

PRIMA C
Article code 6128000

SUPER ALL-ROUNDER

Same meticulously designed optical system as the Prima C but connected to a suspension and swivel arm mounted on a stable base.

The arm can easily be mounted directly on an exam chair.



PRIMA CS
Article code 6128100-801

Labomed Prima CS microscope is also available in two other different versions: 45° head & Ergo head



PRIMA CS, 45° HEAD
Article code 6128100-803



PRIMA CS, ERGO HEAD
Article code 6128100-805

SUPREME CHOICE

Both Prima GN and Prima GN OTS (over the shoulder) are based on a H-shape mobile stand with robust castor wheels. Makes it easy to transport the colposcope to the next room. The unit can be mounted directly on an exam chair.

A powerful LED light source in combination with a superb optical system makes the Prima GN a valuable partner for diagnostic and surgical work.



PRIMA GN, 90° HEAD
Article code 6165000

Labomed Prima GN microscope is also available in two other different versions: 45° head & Ergo head



PRIMA GN, 45° HEAD
Article code 6165000-110



PRIMA GN, ERGO HEAD
Article code 6165000-125

SUPREME CHOICE

Prima GN OTS is also available
in 45° head.



- PRIMA GN OTS, STRAIGHT**
Article code 6165000-100
- PRIMA GN OTS, 45° HEAD**
Article code 6165000-102

LABOMED NUVAR OBJECTIVES

NuVar, adjustable working distance from 200 mm to 270, 300 or 400 mm



NuVar 7
Article code 613340



NuVar 10
Article code 6133401



NuVar 20
Article code 6133351

These objectives provide fluid adjustment of the working distance over several centimeters.

Suitable for Prima OPH and Prima M μ OPH.

Protect NuVar 7, 10 and 20 against splashes, dust and scratches (with a lens protector).



NuVar lens protector 7 & 10
Article code 6156250-801



NuVar lens protector 20
Article code 6133351-812

Double beam Splitter, flat

Article code 6134150



Single Beam Splitter GN

Article code 6122501



If you cannot find what you're looking for, do not hesitate and ask us directly.



Double beam Splitter, flat, kit incl. M200

Article code 6134156




Single Beam Splitter GN version with camera

Article code 6122522



Single Beam Splitter GN version with HDMI camera

Article code 6122525



Double handle for mounting on the magnichanger, GN

Article code 6128600-852

GENERAL SPECIFICATIONS

Product name	PRIMA C	PRIMA CS	PRIMA GN	PRIMA GN OTS
Binocular tube <i>Straight viewing tube 90°</i>	√	√	√	√
Eyepiece <i>WF 10X/18mm with foldable eye guards, diopter adjustment ± 5mm Optional: 12.5X/18mm; fixed eye guards</i>	√	√	√	√
Apochromatic magnichanger <i>5 step > 0.4X, 0.6X, 1.0X, 1.6X and 2.5X</i>	√	√	√	√
Objective <i>f-300, manual fine focus</i>	√	√	√	√
Illumination <i>27W LED, 75,000 Lux</i>	√	√		
Illumination <i>50W LED, 135,000 Lux</i>			√	√
Position LED light source <i>Built in head</i>	√	√		
Position LED light source <i>In swivel arm</i>			√	√
Built in filters <i>Green</i>	√	√		
Built in filters <i>Green & yellow</i>			√	√
Floor mount (*) <i>Star base</i>	√			
Floor mount (*) <i>H-shape</i>		√	√	√

(*) are in the span of delivery, optional a variety is available. Please ask for the options.

All Labomed medical microscopes are treated with MaxLite™. A multi layered coating to guarantee a most efficient transmission of the image. With excellent anti-fog and anti-fungal properties.





Essebaan 50
2908 LK Capelle aan den IJssel
THE NETHERLANDS

+31 (0)10 45 842 22 | info@labomedeuropa.com | www.labomedeuropa.com