


Labomed Lx 300 Gyn/OB a phase contrast microscope for gynecology



A superb microscope to make effortless observations of vaginal smears.

The Lx 300 Gyn/OB is specially composed and comes with a phase contrast set (10x and 40x) and a USB/HDMI camera.

High quality optics in combination with an ergonomic design guarantees effortless observations.

A rechargeable EcoPro LED system delivers day light illumination (recharger is included) Life span of the LED system: approx. 50.000 hours. Total magnification is 1000x with bright field objectives, the phase contrast attachment allows up to 400x magnification.

Article code 9136003-Gyn

