

AXIUS



WORK WITHIN YOUR COMFORT ZONE

LABOMED 



AXIUS

Amazing ergonomics and perfect comfort:

Every dentist's space is designed to provide comfort and efficiency. Medical procedures and patient care require complete harmony between the space, equipment, doctor, assistant and patients. Ergonomic arrangement of devices allows operatories to function smoothly. LABOMED AXIUS enhances safety, comfort, efficiency.

Freedom of movement even in small spaces


AXIUS provides an easy and efficient arrangement of mounting microscope, lights, X-ray equipment, camera monitors and data storage devices all around a central column mounted between the floor and ceiling of the operatory. AXIUS – The axis of efficiency

Center of Comfort


Pain is what an uncomfortable workspace brings. Doctors need to work long hours in arduous procedures and a poorly organized work space causes physical as well as mental fatigue. With all devices within an easy reach of the doctor and positioned ergonomically around the patient, AXIUS is a well-engineered center of command and comfort.

AXIUS

Geometry	Axial four channels to avoid any swivelling collision and safe cable management
Structure	Robust, accident proof metal extrusion with plastic cladding
Diameter	125mm with cladding
Base and ceiling anchoring diameter	230mm
Length availability (for room height)	2.40-2.90 mtr 2.91-3.40 mtr 3.41-3.90 mtr
Weight	10kg / mtr
Load capacity	1250N (125 KgF)
Accessories	Microscope adapter, Monitor arm, Operating light arm & X-ray arm



SPECIFICATIONS



Labomed Europe

Essebaan 50
2908 LK, Capelle aan den IJssel
The Netherlands

Tel: +31 (0) 10 45 842 22

E-mail: sales@labomedeurope.com

www.labomedeurope.com



Distributor

With a policy of continuous development, Labo Europe b.v. reserves the right to change design and specifications without prior notice.

6137120-796/USA/E/1.0 010822