

# ErgoFlex

DENTAL STOOL




# ErgoFlex

Seat Height	Minimum 508mm (20 inches)
Extension	Upto 685mm (27 inches)
Back rest Tilt	0° - 25°
Adjustable Back Rest height	254mm (10 inches)
Arm rest height adjustment	152mm (6 inches)
Arm rest Width Extension	Minimum 533mm (21 inches)
Extension	Upto 736mm (29 inches)
Arm rest radial adjustment	± 50°
Base Diameter	559mm (22 inches)
Fabric Colour	Black
Weight	24 Kg
Seat system design	Contoured water fall design with 110 degree angle
Pneumatic cylinder	Anticorrosive steel, 160 kg capacity
Upholstery	High density PU foam, Anti fading, Anti UV Fabric
Base and arm rest	Made from high polished anti rust Stainless Steel



Armrest Height and Width Adjustment



Armrest Radial & Angular (applies to both armrest)



Backrest Height Adjustment



Seat Height Adjustment



Backrest Tilt Adjustment

## SPECIFICATIONS

### Labomed Europe b.v.

Essebaan 50  
2908 LK Capelle aan den IJssel  
The Netherlands  
info@labomedeuropa.com  
Tel: +31(0)10-4584222  
Fax: +31(0)10-4508251

[www.labomedeuropa.com](http://www.labomedeuropa.com)

With a policy of continuous development, Labo America, Inc. reserves the right to change design and specifications without prior notice.



ISO 9001 : 2008

# LABOMED®

Distributor

6137120-796/USA/E/ 1.0 010821  
© Labo America, Inc. All rights reserved.