

# SLIT LAMP PROTECTIVE SHIELD

SAFETY DEVICE



# PROTECTIVE SHIELD

## YOUR SAFETY IS OUR RESPONSIBILITY

Realizing the present day need for hygiene, LABOMED has designed SPS breath Shield for its Slit lamps to provide relief and safety to Ophthalmologists.



Easy to fit and easy to change, LABOMED safety guard is an ideal solution for providing safety to patient and Doctor at a negligible cost. SPS is built with premium quality plastics that are easily cleanable and can be effectively disinfected. LABOMED SPS provides a well-engineered shield preventing air borne bacteria to travel between the two sides of the instrument. With the use of this breath shield, contamination around the examination table can be greatly reduced.

### FEATURES:

- Easy Installation.
- Easy access to device, Patients & device controls.
- Can fit to any brand Slit Lamp.
- Easy to Disinfect.

### INSTALLATION:

- Remove the Eyepieces.
- Adjust the IPD to Insert the breath Shield.
- Adjust IPD for easy use.
- Peel off the protective film prior to use.



## Labomed Europe b.v.

Essebaan 50  
2908 LK Capelle aan den IJssel  
The Netherlands  
info@labomedeurope.com  
Tel:+31(0)10-4584222  
Fax:+31(0)10-4508251



LABOMED

**LABOMED** 

Distributor

[www.labomedeurope.com](http://www.labomedeurope.com)

With a policy of continuous development, Labo America, Inc. reserves the right to change design and specifications without prior notice.