


# PRIMA ProLine



# THE POWER OF ZOOM



## Assistoscope

For those dental procedures where an assistant is required to be in the same field of view as the doctor, LABOMED Assistoscope with its multiple axis of rotation provides ultimate in comfort and accuracy.




## UViz Fluorescence Lamp

LABOMED Uviz fluorescence lamp attachment is a very convenient attachment for easy identification and differentiation of natural and artificial dental materials. This quick visualisation reduces procedure time and patient discomfort. A UV lamp and barrier filter working in tandem with full safety features provide efficient fluorescence light assisted reparation work on dental carries.



A NEW MILESTONE IN DENTAL MICROSCOPY

# EXTENDERS & ACCESSORIES



Part No.	Description	Part No.	Description
6134115-801	Rotoplate	6134110-800	CCD Attachment
6134150	Double Beam Splitter 50:50		1/2" Camera Mount Adapter
6134160	Double Beam Splitter Inclined 50:50	6134131	DSLR Camera Adapter Kit
6134165	Extender	6134117	Ergo Assistoscope
6214051	Double IRIS Diaphragm	6122015	Ergo Head with Eye Piece
6134120	Adapter Sony (HDR-CX-130) Handycam	6129050	UV Attachment



Labomed Europe  
Essebaan 50  
NL-2908 LK Capelle a/d IJssel  
The Netherlands

tel: +31 (0)10 4584222  
fax: +31 (0)10 4508251  
e-mail: info@labomedeurop.com

[www.labomedeurop.com](http://www.labomedeurop.com)

\*Labomed Europe, a division of Breukhoven b.v.



With a policy of continuous development, Labo America, Inc. reserves the right to change design and specifications without prior notice.

**LABOMED**

Distributor

SPECIFICATIONS

6134150-796 1.1  
© Labo America, Inc. All rights reserved.