# SL 20i/25i





# A PERFECT FIT FOR EYECARE PROVIDERS





A sturdy stand with fluid mechanical movement provides ergonomic benefits for hours of use. Apochromatic optics provide the clearest possible view of the anterior and posterior segments of the patient's eye.

The SL 20i is a nimble system that puts diagnostic essentials at the fingertips of eye care professionals. This microscope delivers all the performance features of the SL 25i in optical and mechanical performance, but with a more compact bottom illumination system without any inclination. This ensures that SL 20i provides doctors with greater space efficiency. The SL 25i delivers ultimate value with a well-designed 2-step optical system for gonioscopy with an efficient 5W LED tower style illumination system inclinable at four steps 5°each.

LABOMED's *MaxLite*<sup>™</sup> coated optics, standard in all SLi series Slit-lamps, yield crisp and detailed images. *MaxLite*<sup>™</sup> coatings provide the highest throughput of light, thanks to superb anti-reflective properties. These optics are anti-fungal, anti-fog and scratch resistant in addition to producing true-to-life resolution.

LABOMED SLi Series Slit lamps are sturdy intruments built with fluidity in movement and ergonomics to benefit hours of use.



The SLi series accessories provide an array of options based on the needs of the doctor's practice. IOP readings are available with the addition of a tonometer. A host of beam splitters and adapters make it possible to integrate third party cameras.



All optics are lead free and all paints used have antimicrobial additives that prevent the growth of bacteria on contact surfaces, ensuring a healthier working environment. All parts in contact with patient are antitoxic and anti-irritation.



"To help doctors for their documentation needs Labomed SL Proline range of camera adapters & smart phone adapters are helpful for still imaging or videography needs"



#### **SL 20i**

Microscope Galilean Magnification 2 step: 10X, 16X

Eyepieces WF 10X/18 mm, dioptre adj. ± 5 mm

12.5x optional

IPD range 49 - 78 mm Slit length 0 - 14 mm

Slit apertures 1, 2, 3, 5,8,10 & 14 mm

Slit rotation 0 - 180° Slit tilt

Red free & Cobalt blue **Filters** 

Magnification Obj. Eye pc Mag FOV(mm) 1 12.5 10 22.5 1.6 12.5 16 14.1 10 8 22.5 1 1.6 10 12.8 14.1

Movement ranges:

Longitudinal (in/out) 80 mm Lateral (left/right) 109 mm Vertical (up/down) 30 mm Chin rest range 80 mm Illumination 5W LED

367 mm X 360 mm X 5.5, Metal Table base size

Weight 25 kg (packed)

Optional accessories Tonometer, Breath shield, Beam splitter, Camera adapter, Diffuser, Motorized table, Smart phone adapters

	SL 20i	SL 25i
All dimensions in mm		

#### **DESCRIPTION**

SL 20i LED SL 25i LED

#### SL25i

Galilean 2 step: 10X, 16X

WF 10X/18 mm, dioptre adj. ± 5 mm

12.5x optional 49 - 78 mm 0 - 12 mm

0.3, 1, 3, 5, 9 & 12 mm

0 - 180°

20° (Four inclination steps 5° each) Red free, Cobalt blue, Heat absorbing,

Neutral density

Obj. Eye pc Mag FOV(mm) 1 12.5 10 22.5 1.6 12.5 16 14.1 8 22.5 1 10 1.6 10 12.8 14.1

80 mm 109 mm 30 mm 80 mm 5W LED

367 mm X 360 mm X 5.5, Metal

30 kg (packed)

#### CONFIGURATION CATALOG NO.

2 step Slit Lamp 8126701 2 step Slit Lamp 8127701



### Labomed Europe b.v.

Essebaan 50 2908 LK Capelle aan den IJssel The Netherlands info@labomedeurope.com Tel:+31(0)10-4584222 Fax:+31(0)10-4508251



## www.labomedeurope.com





Distributor

© Labo America, Inc. All rights reserved 8126701-796 1.0 011019