




Datasheet

Sterile microscope drapes CMO Ø65 mm, pack of 20 pcs

Article code: A612099-527

Product Details

- **Size (HxW):** 117 x 163 cm
- **Quantity:** 20
- **Contains latex:** No
- **Sterile:** Yes
- **Sterilization method:** EO
- **Material:** Polyethylene
- **Intended use:** Equipment drapes are used as a protective covering for various surgical and/or nonsurgical instruments/equipment [e.g., a surgical or endoscopic tray cover or liner, an operating room (OR) equipment]
- **Single use:** Yes



A spacer ring is required to compensate for the diameter difference between the drape (Ø65 mm) and the microscope CMO (Ø59 mm).

Article code: A612099-529

Regulatory & Compliance

- **CE Notifying Body:** CE0044
- **Medical Classification:** IS
- **GMDN Code:** 12535
- **DEHP/Phthalates:** No
- **Bisphenol A:** No
- **PVC:** No
- **REACH Compliance:** Yes
- **IFU:** Yes

Packaging & Shelf Life

- **Box Size (LxWxH mm):** 455 x 295 x 290
- **Box M³:** 0.0389
- **Box Weight (g):** 4450
- **Cases Per Pallet Conversion:** 39
- **Shelf Life:** 5 Years
- **Double Sterilization Allowed:** No

Information Last Updated: 21 Jun 2022