



LABOMED[®]
ideas for vision

PRIMAlite

Labomed dental microscopes meet the most demanding expectations because they combine in the same equipment everything you expect:

- Big in features
- Big in finishings
- Big in adaptability
- Big in customer care
- Big in competitive pricing, all included and without surprises.

Think big



Think big

PRIMAlite

At Labomed we know that investing in a dental microscope for your practice is an important step towards the quality and service you offer to your patients and also towards the excellence and prestige of your center. For this reason, it is a decision that must be made with perspective, since it is an asset that you will have for years. So, think big.

Labomed offers the perfect dental microscope. Exactly the one responding to your demands and expectations. Everything you desire combined in one microscope: great performance, great finishes, great adaptability, great customer service and a great price, all included and without surprises.

In short, with Labomed solutions you will have a high-end dental microscope, adapted to your specific needs and with a team of people who will be at your side at all times.

Think big. What features do you need, what type of lenses? Do you want to capture images or project the image on a screen? Should the microscope be mounted on the ceiling or wall or should it be on wheels? Think big. With Labomed it is possible.



PRIMA^{lite}

Think big



Big in features

Labomed **PRIMA lite** dental microscopes come equipped as standard with unique, high-tech features designed by and for dentists. Apochromatic optics. Unique, intense and shadowless lighting. Lens protector. Ergonomics. Smooth control and handling.



Big in finishings

At Labomed, quality is our obsession. Our microscopes are made with the highest quality materials. The equipment comes with antibacterial paint and a MaxLite coating on the optics, prisms and mirrors. And, in addition, they offer a modern and impeccable look, just like your cabinet deserves.



Big in adaptability

Labomed offers a choice of multiple configurations. Always with your needs in mind. Choose the model. Determine whether it should be a fixed installation, on the ceiling or on the wall or on the Axis column. Or should it be a mobile stand.

Do you wish to capture images with a camera or with a smartphone? More than enough options.

For us, your wish always comes first.



Big in customer care

At Labomed, collaboration is of paramount importance.

Our involvement is present throughout the sales process and we provide support where required. Human contact always comes first. Pre-sales and after-sales, from start to finish!



Big in competitive pricing

Nobody like Labomed offers you so much in a single price without surprises or unnecessary extras: great features, great finishes, great adaptability and great customer service.

PRIMA^{lite}



PRIMA^{lite}

Specifications

Observation tubes	0-210° inclinable binocular
Type of eye guards	Foldable rubber
Eyepieces standard	WF 10x/18 mm
Option	WF 12.5x/18 mm
Diopter adjustment	± 7
IPD	49-78 mm
Magnichanger	5 step
Magnichanger steps	0.4x, 0.6x, 1.0x, 1.6x and 2.5x
Objective standard	f 250 mm
Objective optional	f 300 mm, f 400 mm, NuVarTM 7 or NuVarTM 10
Light source	27 W
LED light	Yes
Lux power	100.000
Light source location	Head
Optical fiber	No
Minimal LED lifetime	50.000 hours
Built-in filters	Green and yellow
Fluent carrier	Yes

And additionally, all **PRIMA lite** microscopes come equipped with:

- Antibacterial paint.
- MaxLite coating on optics, prisms and mirrors.
- Apochromatic optics.
- Lens protector that keeps the optics clean from splashes of spraying and other waste.
- Automatic light off.
- 5+ microscope mounting options and interchangeable.
- Set of covers for the controls, suitable for the autoclave and respectful with the environment by avoiding creating more waste.



Let's think big together!

Contact us here:

Eugene C. Dams
Sales Manager North, West & East Europe
Labomed Europe b.v.
+31 624693712
eugene@labomedeuropa.com

Eugenia Clavero
Sales Manager South Europe
+31 683586788
eugenia@labomedeuropa.com

Labomed Europe
T +31(0)10 45 842 22
Essebaan 50
2908 LK Capelle aan den IJssel
The Netherlands
www.labomedeuropa.com